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# English

## Measurement of the migration of low-molecular substances from packaging's by sample preparation with the migration cell System Siegwerk (Sieg-Mi-Flex) and following analytical quantification.

In the EU, for packaging with direct food contact, the standards of the EU regulation 1935/2004, Article 3, apply. The most important requirement basically states that all materials and objects shall be manufactured as to assure that their material components will only migrate in such small quantities that any possible health risks can be excluded.

The migration cell System Siegwerk (Sieg-Mi-Flex) simplifies sample preparation while dealing with migration analytics. With these migration cells, it is possible to achieve a migration estimation of low-molecular substances, for example of food packaging, that takes place even right at the ready-made rolls. Owing to the flexibility of cell system, all sorts of samples can be analyzed with the correct choice of the liquid food simulant. After the preparation and storage of the cells, according to the EU regulations, the simulant solutions can be quantified via the respective analytic methods (GC/MS, HPLC-MS, HPLC-DAD oder DC), or after a prior concentration of samples via SPE, SPME or rotary evaporation.

The migration cell System Siegwerk (Sieg-Mi-Flex) includes a ring shaped volume element (center ring) that has sealing elements that are arranged at it's frontal area (O-rings made of silicone, surrounded by FEP). The circular area is interrupted by two openings, that are filled with one filler pipe each. The filler pipe is eather sealed with a common laboratory screw with a cap and a GL thread or with a septum. The volume element is being braced between the front and the bottom plate by four screws, with the bolts having ergonomically formed handwheels. The circular volume element (center ring) is offered in different sizes, making it possible that the migration cell can be adjusted rapidly to the current trial area.

Because of the modular construction of the migration cell (Sieg-Mi-Flex), an adaption of the cell for fixed adsorbents can be expected in the future.


## Stainless steel migration cell System Siegwerk (Sieg-Mi-Flex)

A transition of low-molecular substances from packaging to the filling material can take place in two different ways:

1. In the batch or in the roll by set off.
2. The migration of substances from the printed ink side through the foil material into the filling goods.

For migration analysis of various samples, Eva Holster, Dr. Dieter Franke (Siegwerk Druckfarben AG & Co. KGaA) and Werner Zillger (LABC Labortechnik, Hennef) developed the System Siegwerk(Sieg-Mi-Flex). It guarantees the highest possible flexibility regarding different surface sizes and structures of the samples.

1. Central rings with smaller diameters than nom. width DN120 achieve flexibility in analysis of sample sizes. A fixing plate clamps the central ring into the housing made of stainless steel.

2. Ergonomically formed plastic hand-wheels ensure an easy handling in locking and opening of the stainless steel migration cells. Additionally, the compact and space saving design of these cells allows horizontal or vertical storage e.g. for tempering in an oven.

3. Flexible, chemical resistant FEP/Silicon O-rings are used for sealing elements. The diameter of the filler necks is big enough for a practicable handling in filling the volume of the central ring. The volume ring can be closed/locked by temperature stable laboratory standard GL14 cylindrical cup.

Eva Holster and Dr. Dieter Franke, member of analytical department Siegwerk Druckfarben AG & Co. KGaA, located in Siegburg remarked: „The migration cell “System Siegwerk” offers an optimum of sealing between sample and cells, a maximum convinienis in daily handling and guarantees the highest possible accuracy in migration analysis“.

### Technical data:

The migration cell “**System Siegwerk**” consists of:

- One Fixing plate (A) and different central rings (B) of your choice.

#### A) Fixing plate (fits all following central rings!):

- Stainless steel (VA 1.4571) plates with threaded holes, inside polishes,
- 4 x stainless steel bolt, plastic star nut and thread (DIN 125B A2)

**Part No.: 715318: Price on request!**

#### B) Central rings:

- Stainless steel (VA 1.4571) in and filling connecting pieces with GL14 thread inclusive. Cover screw-cap GL14 (PBT housing + PTFE/Sil. - Seal)
- 2 x O-ring FEP/Silicon suitable for groove in central ring (top and bottom)
- Aluminium-positioning-plate (except for the nom. width DN120!)

Central rings article No.	Nom. width	One-sided approx. contact area	Two-sided approx. contact area	Capacity (ml)
<b>715319</b>	DN120	1.00 dcm <sup>2</sup>	2.0 dcm <sup>2</sup>	200
<b>715320</b>	DN100	0.75 dcm <sup>2</sup>	1.5 dcm <sup>2</sup>	150
<b>715321</b>	DN80	0.50 dcm <sup>2</sup>	1.0 dcm <sup>2</sup>	100
<b>715322</b>	DN70	0.35 dcm <sup>2</sup>	0.7 dcm <sup>2</sup>	70
<b>715323</b>	DN60	0.25 dcm <sup>2</sup>	0.5 dcm <sup>2</sup>	50
<b>715324</b>	DN30	0.05 dcm <sup>2</sup>	0.1 dcm <sup>2</sup>	10

**Prices on request!**


## Glass-Migration-Cell System-Siegwerk (Sieg-Mi-Flex)

used for sample preparation in the measurement of migration of low molecular weight substances from packaging materials and subsequent analytical quantification.

A transition of low-molecular substances from packaging to the filling material can take place in two different ways:

1. In the batch or in the roll by set off.
2. The migration of substances from the printed ink side through the foil material into the filling goods.

For migration analysis of various samples, Eva Holster, Dr. Dieter Franke (Siegwerk Druckfarben AG & Co. KGaA) and Werner Zillger (LABC Labortechnik, Hennef) developed the glass-migration-cell "System Siegwerk".

1. Middle rings from the same high made from borosilicate glass 3.3, with the nominal diameter of DN120 with the help of two plan glass-fixing-plates clamped.
2. Ergonomically formed plastic hand-wheels ensure an easy handling in locking and opening of the glass steel migration cell. Additionally, the compact and space saving design of these cells allow horizontal or vertical storage e.g. for tempering in an oven.
3. Flexible, chemical resistant FEP/Silicon O-rings are used for sealing elements. The diameter of the filler necks is big enough for a practicable handling in filling the volume of the central ring. The volume ring can be closed/locked by temperature stable laboratory standard GL14 cylindrical cup.

### Technical data from the *glass-migration-cell System-Siegwerk* :

Two fixing plates and one central ring with Nom. width (DN120), One-sided approx. contact area = 1,00 dm<sup>2</sup>, Two-sided approx. contact area = 2,0 dm<sup>2</sup>, Capacity 200ml.

- 2 borofloat plates with holes, protective inserts for screws, bolts and nuts
- 1 borosilicate glass ring with O-ring DN120-slots, 2 x GL14, ID max. 8.5 to 9 mm,

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