

Siebgewebe aus ETFE

Artikelnummer	Länge des Siebgewebes	Maschenweite	offene Siebfläche	Gewicht	Dicke	
717946-1	1m	70µm	21,75%	123g/cm ²	158µm	
717946-3	3m					
717946-10	10m					
717947-1	1m	85µm	27%	125g/cm ²	175µm	
717947-3	3m					
717947-10	10m					
717948-1	1m	105µm	32,3%	103g/cm ²	155µm	
717948-3	3m					
717948-10	10m					
717949-1	1m	120µm	36%	113g/cm ²	170µm	
717949-3	3m					
717949-10	10m					
717950-1	1m	150µm		110g/cm ²	200µm	
717950-3	3m					
717950-10	10m					
717951-1	1m	210µm	46%	88g/cm ²	172µm	
717951-3	3m					
717951-10	10m					
717952-1	1m	250µm	39%	165g/cm ²	290µm	
717952-3	3m					
717952-10	10m					
717953-1	1m	300µm	36%	230g/cm ²	385µm	
717953-3	3m					
717953-10	10m					
717954-1	1m	420µm	39%	260g/cm ²	450µm	
717954-3	3m					
717954-10	10m					
717955-1	1m	500µm	42%	270g/cm ²	500µm	
717955-3	3m					
717955-10	10m					
717956-1	1m	590µm	47%	240g/cm ²	485µm	
717956-3	3m					
717956-10	10m					
717957-1	1m	850µm	57,6%	180g/cm ²	505µm	
717957-3	3m					
717957-10	10m					
717958-1	1m	1000µm	44,5%	480g/cm ²	1000µm	
717958-3	3m					
717958-10	10m					