

Die neue Summit-Reihe

Modell	Ablesbar-keit	Wägebereich	Waagschale	Reproduzierbarkeit	Meßzeit	Bauform		Bestellnummer
Summit-S (mit externem Kalibriergewicht)								
S-234	0,1 mg	230 g	∅ 90mm	0,1 mg	2,5 sec	1		717980
S-114	0,1 mg	110 g	∅ 90 mm	0,1 mg	2,5 sec	1		717981
S-64	0,1 mg	60 g	∅ 90 mm	0,1 mg	2,5 sec	1		717982
S-603	0,001 g	600 g	∅115 mm	0,001 g	1 sec	2		717983
S-403	0,001 g	400 g	∅115 mm	0,001 g	1 sec	2		717984
S-203	0,001 g	200 g	∅ 115 mm	0,001 g	1 sec	2		717985
S-6002	0,01 g	6000 g	180 x 180 mm	0,01 g	1 sec	3		717986
S-4002	0,01 g	4000 g	180 x 180 mm	0,01 g	1 sec	3		717987
S-2002	0,01 g	2000 g	180 x 180 mm	0,01 g	1 sec	3		717988
S-602	0,01 g	600 g	∅ 115 mm	0,01 g	1 sec	4		717989
S-402	0,01 g	400 g	∅ 115 mm	0,01 g	1 sec	4		717990
S-8001	0,1 g	8000 g	180 x 180 mm	0,1 g	1 sec	3		717991
S-4001	0,1 g	4000 g	180 x180 mm	0,1 g	1 sec	3		717992

Justiergewichte

Modell	Genauigkeitsklasse	Gewicht	Bestell-Nr.	
S-234	E2	200 g	870200.8	
S-114	E2	100 g	870100.8	
S-64	E2	50 g	870050.8	
S-603	E2	500 g	870500.8	
S-403	E2	200 g	870200.8	
S-203	F1	200 g	870200.5	
S-6002	E2	5000 g	875000.8	
S-4002	E2	2000 g	872000.8	
S-2002	F2	2000 g	872000.6	
S-602	F2	600 g	870500.6	
S-402	F2	400 g	870200.6	
S-8001	F2	5000 g	875000.6	
S-4001	F2	2000 g	872000.6	