

LABC

LABORTECHNIK



LABC

MOSH/MOAH - KIT

Number: 120-0962797

Analytik mit ZeroSept®AIR

The ZeroSept®AIR septum from LABC-Labortechnik provides a leak-proof and mineral oil-free solution to avoid contamination during automated MOSH/MOAH analysis. The LABC kit connects the septum to a fine thread vial to perform MOSH/MOAH analyses in an interference-free environment.

More info:

Phone: 02242-969460

Email: info@labc.deWebsite: www.labc.de

Parts

Vial



Fine thread vial

The 10ml fine thread vial is made of clear glass of the first hydrolytic class. It also has a rounded base.

Art. Nr.: 120-0588751

Kappe



Fine thread screw cap

A magnetic 18mm ZeroSept®AIR closure consisting of steel and a 2.5mm thick septum made of PTFE.

Art. No.: 120-0916660