

**LABC**

LABORTECHNIK



# LABC-

## MICROLITER VIAL BOTTLES

### Microliter Vials

Our microliter bottles are the absolute specialty vials for precious and small sample quantities.

They are precisely manufactured to ensure optimum residual emptying. The microliter vials are "Made in Germany".

#### More infos:



Phone: +49-2242-96946-0

Email: [info@labc.de](mailto:info@labc.de)

Website: [www.labc.de](http://www.labc.de)

### Our variants

**iV2 $\mu$** 

#### **iV2 $\mu$ -Vial - Insert - Vial**

Slim drawn and directly closable insert vial with a volume of 250 $\mu$ l.

**Vase**

#### **Vase-Vial**

The vase vial has a volume of up to 1.2ml and a tapered base for optimal residual emptying.

**i3 $\mu$** 

#### **i3 $\mu$ -Vial $\mu$ -Vial**

The  $\mu$ -Vial with integrated 350  $\mu$ l insert. The fused insert provides a common bottom with the outer vial.